


Anti-Mmp2 Rabbit Monoclonal Antibody For Wb And Elisa - P08253

Item Number: CM0006484



Introduction

Rabbit monoclonal antibody against human MMP2 (72 kDa type IV collagenase, gelatinase A, CLG4), validated for WB and ELISA. Detects human, mouse, rat MMP2. Research use in angiogenesis, tumor invasion, and matrix remodeling.

[Learn More](#)

Parameter	Value
Product Name	MMP2 Rabbit mAb
Remarks/Alias	CLG4; MMP2
Species	Human
Gene ID	4313
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 451-660 of human MMP2 (NP_004521.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08253
Protein Weight	74 kDa
Shipping	Ice bag