


Anti-Mmp3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08254

Item Number: CM0018145



Introduction

MMP3 Rabbit pAb validated for WB, IHC-P, IF/ICC, ELISA in Human, Mouse, Rat. Targets stromelysin-1, a key matrix metalloproteinase in ECM remodeling, inflammation, and antiviral immunity. Ideal for cancer and arthritis research.

[Learn More](#)

Parameter	Value
Product Name	MMP3 Rabbit pAb
Remarks/Alias	SL-1; STMY; STR1; CHDS6; MMP-3; STMY1; MMP3
Species	Human
GeneID (Human)	4314
GeneID	4314
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08254
Protein Weight	54kDa
Shipping	Ice bag