

Anti-Mnx1/Hb9/Hlxb9 Polyclonal Antibody For Wb And Elisa - P50219

Item Number: CM0021424

Introduction

Rabbit polyclonal antibody targeting human MNX1/HB9/HLXB9 transcription factor. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for motor neuron and pancreas development studies.

[Learn More](#)

Parameter	Value
Product Name	MNX1/HB9/HLXB9 Rabbit pAb
Remarks/Alias	HB9; HLXB9; SCRA1; HOXHB9; MNX1/HB9/HLXB9
Species	Human
GeneID (Human)	3110
GeneID	3110
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 201-315 of human MNX1/HB9/HLXB9 (NP_005506.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50219
Protein Weight	41kDa
Shipping	Ice bag