

Anti-Mpi Rabbit Polyclonal Antibody For Wb, Elisa - P34949

Item Number: CM0026128

Introduction

High-quality rabbit polyclonal antibody against human MPI, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying glycosylation and congenital disorders of glycosylation.

[Learn More](#)

Parameter	Value
Product Name	MPI Rabbit pAb
Remarks/Alias	PMI; PMI1; CDG1B; MPI
Species	Human
GeneID (Human)	4351
GeneID	4351
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-423 of human MPI (NP_002426.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34949
Protein Weight	47kDa
Shipping	Ice bag