


Anti-Mrp3 Rabbit Polyclonal Antibody For If-P, Elisa - O15438

Item Number: CM0019020



Introduction

High-quality anti-MRP3 rabbit polyclonal antibody, validated for IF-P and ELISA. Detects human MRP3 (ABCC3, 169 kDa), an ABC transporter involved in multidrug resistance, bile acid transport, and hepatocellular export. Ideal for liver and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MRP3 Rabbit pAb
Remarks/Alias	MLP2, MRP3, ABC31, MOAT-D, cMOAT2, EST90757, ABCC3
Species	Human
Gene ID	8714
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 855-954 of human MRP3 (NP_003777.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human
SWISS	O15438
Protein Weight	169kDa
Shipping	Ice bag