


Anti-Mrpl12 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P52815

Item Number: CM0005312



Introduction

High-specificity anti-MRPL12 rabbit monoclonal antibody validated for WB, IF/ICC, IHC-P, and ELISA. Detects human, mouse, and rat MRPL12, a mitochondrial ribosomal protein involved in translation and transcription. Ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	MRPL12 Rabbit mAb
Remarks/Alias	5c5-2; L12mt; MRPL7; RPML12; MRPL7/L12; MRP-L31/34
Species	Human
GeneID	6182 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 46-198 of human MRPL12 (NP_002940.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52815
Protein Weight	21 kDa
Shipping	Ice bag