


## Anti-Mrpl20 Polyclonal Antibody For Wb And Elisa - Q9Byc9

Item Number: CM0021619



### Introduction

Anti-MRPL20 polyclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat. Ideal for studying mitochondrial ribosome structure and protein synthesis. SWISS: Q9BYC9.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] MRPL20 Rabbit pAb                                                                                          |
| Remarks/Alias    | L20mt; bL20m; MRP-L20; MRPL20                                                                                             |
| Species          | Human                                                                                                                     |
| Gene ID          | 55052 (Human)                                                                                                             |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-149 of human MRPL20 (NP_060441.2) |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, ELISA                                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | Q9BYC9                                                                                                                    |
| Protein Weight   | 17kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |