


Anti-Mrpl3 Polyclonal Antibody For Wb, Elisa - P09001

Item Number: CM0024175



Introduction

Rabbit polyclonal antibody targeting human mitochondrial ribosomal protein L3 (MRPL3). Suitable for Western blotting and ELISA; validated in human samples. Ideal for mitochondrial ribosome and protein synthesis research.

[Learn More](#)

Parameter	Value
Product Name	MRPL3 Rabbit pAb
Remarks/Alias	MRPL3; COXPD9; MRL3; RPML3; mitochondrial ribosomal protein L3
Species	Human
Gene ID (Human)	11222
Gene ID	11222
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 249-348 of human MRPL3 (NP_009139.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09001
Protein Weight	38kDa
Shipping	Ice bag