

## Anti-Mrpl36 Rabbit Polyclonal Antibody For Wb And Elisa - Q9P0J6

Item Number: CM0023535



### Introduction

High-quality rabbit polyclonal antibody against human MRPL36 (bL36m), validated for Western blot and ELISA. Ideal for mitochondrial ribosome research and protein synthesis studies.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | MRPL36 Rabbit pAb                                                                                                  |
| Remarks/Alias    | RPMJ, BRIP1, L36mt, bL36m, PRPL36, MRP-L36                                                                         |
| Species          | Human                                                                                                              |
| GeneID           | 64979                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 39-103 of human MRPL36 (NP_115868.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human                                                                                                              |
| SWISS            | Q9P0J6                                                                                                             |
| Protein Weight   | 12kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |