


Anti-Mrpl52 (Kd Validated) Rabbit Polyclonal Antibody For Wb And Elisa - Q86Ts9

Item Number: CM0020541



Introduction

KD Validated anti-MRPL52 rabbit polyclonal antibody for Western blot and ELISA. Detects human mitochondrial ribosomal protein L52 (14 kDa, Q86TS9). Suitable for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRPL52 Rabbit pAb
Remarks/Alias	mL52, 39S ribosomal protein L52, L52mt, MRP-L52, Mrpl52, RM52
Species	Human
GeneID (Human)	122704
GeneID	122704
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 24-123 of human MRPL52 (NP_848026.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86TS9
Protein Weight	14kDa
Shipping	Ice bag