


Anti-Mrps25 Rabbit Polyclonal Antibody For Wb, Elisa - P82663

Item Number: CM0027015



Introduction

CamelBio offers Rabbit Polyclonal Anti-MRPS25 Antibody for WB and ELISA on human samples. Targets 20kDa mitochondrial ribosomal protein S25, relevant to mitochondrial disease research and assay development.

[Learn More](#)

Parameter	Value
Product Name	MRPS25 Rabbit pAb
Remarks/Alias	RPMS25; COXPD50; MRP-S25; MRPS25
Species	Human
GeneID (Human)	64432
GeneID	64432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 104-173 of human MRPS25 (NP_071942.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P82663
Protein Weight	20kDa
Shipping	Ice bag