


Anti-Mtap Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q13126

Item Number: CM0002668



Introduction

Rabbit monoclonal antibody targeting human MTAP (S-methyl-5'-thioadenosine phosphorylase). Validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse, rat. Suitable for polyamine metabolism and cancer research. MW 31kDa.

[Learn More](#)

Parameter	Value
Product Name	MTAP Rabbit mAb
Remarks/Alias	BDMF; MSAP; DMSFH; LGMBF; DMSMFH; c86fus; HEL-249; MTAP
Species	Human
Gene ID (Human)	4507
Gene ID	4507
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 184-283 of human MTAP (Q13126).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13126
Protein Weight	31kDa
Shipping	Ice bag