

Anti-Mtfmt Rabbit Polyclonal Antibody For Wb, Elisa - Q96Dp5

Item Number: CM0021893

Introduction

Rabbit polyclonal antibody against human MTFMT (mitochondrial methionyl-tRNA formyltransferase). Recognizes human, mouse, and rat. Validated for Western blot and ELISA. Ideal for mitochondrial translation and disease research.

[Learn More](#)

Specification	Detail
Product Name	MTFMT Rabbit pAb
Remarks/Alias	FMT1; COXPD15; MC1DN27; MTFMT
Species	Human
GeneID	123263
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-389 of human MTFMT (NP_640335.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DP5
Protein Weight	44kDa
Shipping	Ice bag