

Anti-Mthfd1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P11586

Item Number: CM0027049



Introduction

MTHFD1 rabbit polyclonal antibody for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat. Targets cytoplasmic C-1-tetrahydrofolate synthase, a trifunctional enzyme in one-carbon metabolism and nucleotide biosynthesis.

[Learn More](#)

Parameter	Specification
Product Name	MTHFD1 Rabbit pAb
Remarks/Alias	CIMAH; MTHFC; MTHFD; MTHFD1
Species	Human
GeneID (Human)	4522
GeneID	4522
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human MTHFD1 (NP_005947.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11586
Protein Weight	102kDa
Shipping	Ice bag