


Anti-Mtmr8 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Ef0

Item Number: CM0023500



Introduction

MTMR8 rabbit polyclonal antibody validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, and rat. Targets MTMR8 (Q96EF0), a lipid phosphatase regulating autophagy. Ideal for autophagy and lipid signaling research.

[Learn More](#)

Parameter	Value
Product Name	MTMR8 Rabbit pAb
Remarks/Alias	MTMR8
Species	Human
Gene ID (Human)	55613
Gene ID	55613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 430-560 of human MTMR8 (NP_060147.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EF0
Protein Weight	79kDa
Shipping	Ice bag