

Anti-Mtx2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O75431

Item Number: CM0007338

Introduction

Rabbit monoclonal antibody targeting human MTX2, a mitochondrial import protein. Suitable for WB, IHC-P, ELISA. Reacts with human, mouse, rat; detects ~30 kDa band. Ideal for mitochondrial function studies and IVD development.

[Learn More](#)

Parameter	Value
Product Name	MTX2 Rabbit mAb
Remarks/Alias	MDPS; metaxin-2
Species	Human
GeneID (Human)	10651
GeneID	10651
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-263 of Human MTX2 (NP_006545.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75431
Protein Weight	30 kDa
Shipping	Ice bag