

Anti-Muc1 Rabbit Polyclonal Antibody For Wb, Elisa - P15941

Item Number: CM0023174



Introduction

High-specificity anti-MUC1 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat MUC1. Ideal for cancer biomarker research and epithelial cell studies.

[Learn More](#)

Parameter	Value
Product Name	MUC1 Rabbit pAb
Remarks/Alias	EMA; MCD; PEM; PUM; KL-6; MAM6; MCKD; PEMT; CD227; H23AG; MCKD1; MUC-1; ADMCKD; ADTKD2; Ca15-3; ADMCKD1; CA 15-3; MUC-1/X; MUC1/ZD; MUC-1/SEC; MUC1
Species	Human
GeneID (Human)	4582
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human MUC1 (P15941).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15941
Protein Weight	122kDa
Shipping	Ice bag