

Anti-Muc20 Rabbit Polyclonal Antibody For Wb And Elisa - Q8N307

Item Number: CM0023450

Introduction

High-quality Rabbit Polyclonal Antibody against human MUC20, validated for Western blot (WB) and ELISA. Recognizes human, mouse, rat MUC20, a regulator of MET signaling. Ideal for kidney and urogenital cancer research.

[Learn More](#)

Parameter	Value
Product Name	MUC20 Rabbit pAb
Remarks/Alias	MUC-20; MUC20
Species	Human
GeneID (Human)	200958
GeneID	200958
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 100-200 of human MUC20 (NP_689886.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N307
Protein Weight	72kDa
Shipping	Ice bag