

Anti-Mxi1 Rabbit Polyclonal Antibody For Wb, Elisa - P50539

Item Number: CM0027859

Introduction

Rabbit polyclonal antibody against human MXI1, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying MYC transcriptional antagonism and nucleosomal regulation.

[Learn More](#)

Parameter	Value
Product Name	MXI1 Rabbit pAb
Remarks/Alias	MXI; MAD2; MXD2; bHLHc11; MXI1
Species	Human
GeneID (Human)	4601
GeneID	4601
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 76-295 of human MXI1 (NP_569157.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50539
Protein Weight	26kDa
Shipping	Ice bag