

Anti-Myd88 Monoclonal Antibody For Wb, Elisa - Q99836

Item Number: CM0006860



Introduction

Rabbit monoclonal antibody against human MyD88, validated for WB and ELISA with cross-reactivity to mouse and rat. Essential for studying TLR/IL-1R signaling and innate immunity research applications.

[Learn More](#)

Parameter	Value
Product Name	MyD88 Rabbit mAb
Remarks/Alias	WM1; IMD68; MYD88D; MyD88
Species	Human
Gene ID (Human)	4615
Gene ID	4615
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human MyD88 (NP_002459.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99836
Protein Weight	33kDa
Shipping	Ice bag