

Anti-Myeloperoxidase (Mpo) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P11247

Item Number: CM0022233

Introduction

Rabbit polyclonal antibody targeting mouse myeloperoxidase (Mpo), suitable for WB, IF-P, IHC-P, and ELISA. Detects Mpo in human, mouse, and rat samples, enabling research on neutrophil biology, NETosis, and oxidative immune responses.

[Learn More](#)

Parameter	Value
Product Name	Mouse Mpo Rabbit pAb
Remarks/Alias	mKIAA4033; Mouse Mpo
Species	Mouse
GeneID (Human)	4353
GeneID	17523
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-718 of mouse Mpo (NP_034954.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11247
Protein Weight	81kDa
Shipping	Ice bag