

Anti-Myl12B Monoclonal Antibody For Wb, If/Icc, Elisa - O14950

Item Number: CM0009058

Introduction

Recombinant rabbit monoclonal antibody against MYL12B (MLC-2A, MRLC2). Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for cytoskeleton and cell motility research.

[Learn More](#)

Parameter	Value
Product Name	MYL12B Rabbit mAb
Remarks/Alias	MLC-B; MRLC2; MYL12B
Species	Human
GeneID (Human)	103910
GeneID	103910
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 73-172 of human MYL12B (O14950).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14950
Protein Weight	20kDa
Shipping	Ice bag