

Anti-Myl2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P10916

Item Number: CM0015303



Introduction

MYL2 Rabbit pAb specifically detects ventricular/cardiac myosin regulatory light chain 2 in human, mouse, and rat. Ideal for WB, IF-P, and ELISA applications in cardiac muscle contraction and cardiomyopathy research.

[Learn More](#)

Parameter	Value
Product Name	MYL2 Rabbit pAb
Remarks/Alias	MLC2; CMH10; MFM12; MLC-2; MLC-2v; MLC-2s/v; MYL2
Species	Human
GeneID (Human)	4633
GeneID	4633
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-166 of human MYL2 (NP_000423.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P10916
Protein Weight	19kDa
Shipping	Ice bag