

Anti-Myo1D Polyclonal Antibody For Wb And Elisa - 094832

Item Number: CM0020259

Introduction

Rabbit polyclonal antibody against human MYO1D for WB and ELISA. Detects human and mouse unconventional myosin-I δ (116 kDa), a motor protein involved in endosomal trafficking and ciliary polarity.

[Learn More](#)

Parameter	Description
Product Name	MYO1D Rabbit pAb
Remarks/Alias	myr4; PPP1R108; MYO1D
Species	Human
GeneID (human)	4642
GeneID	4642
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-380 of human MYO1D (NP_001290209.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	094832
Protein Weight	116kDa
Shipping	Ice bag