

Anti-N-Myc/Mycn Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P04198

Item Number: CM0013299



Introduction

Validate N-Myc/MYCN expression in neuroblastoma and neuronal cells with this rabbit polyclonal antibody, suitable for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat species.

[Learn More](#)

Parameter	Value
Product Name	N-Myc/MYCN Rabbit pAb
Remarks/Alias	NMYC; ODED; MODED; N-myc; bHLHe37; MYCNsORF; MYCNsPEP; N-Myc/MYCN
Species	Human
Gene ID (Human)	4613
Gene ID	4613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 168-267 of human N-Myc/MYCN (NP_005369.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04198
Protein Weight	50kDa
Shipping	Ice bag