


Anti-Naa25 Polyclonal Antibody For Wb And Elisa - Q14Cx7

Item Number: CM0029941



Introduction

Rabbit polyclonal antibody against NAA25 (NatB auxiliary subunit). Suitable for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying protein N-terminal acetylation and cell-cycle progression.

[Learn More](#)

Parameter	Value
Product Name	NAA25 Rabbit pAb
Remarks/Alias	NAP1; MDM20; C12orf30; NAA25
Species	Human
GeneID (Human)	80018
GeneID	80018
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 703-972 of human NAA25 (NP_079229.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14CX7
Protein Weight	112kDa
Shipping	Ice bag