


Anti-Naa40 Polyclonal Antibody For Wb, Elisa - Q86Uy6

Item Number: CM0028721



Introduction

High-quality rabbit polyclonal antibody specific to human NAA40 (NatD), an N-alpha-acetyltransferase for histones H4 and H2A. Validated for Western blot and ELISA. Research use only.

[Learn More](#)

Parameter	Value
Product Name	NAA40 Rabbit pAb
Remarks/Alias	NatD; NAT11; PATT1; hNatD; NAA40
Species	Human
Gene ID (Human)	79829
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-237 of human NAA40 (NP_079047.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86UY6
Protein Weight	27kDa
Shipping	Ice bag