


Anti-Nacc1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Re7

Item Number: CM0015356



Introduction

Anti-NACC1 rabbit polyclonal antibody for WB and ELISA applications. Cross-reacts with human and mouse NACC1. Suitable for research on transcriptional repression and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	NACC1 Rabbit pAb
Remarks/Alias	NAC1; BEND8; NAC-1; NECFM; BTBD30; BTBD14B; NACC1
Species	Human
GeneID (Human)	112939
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 170-340 of human NACC1 (NP_443108.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96RE7
Protein Weight	57kDa
Shipping	Ice bag