

Anti-Nadph Oxidase 4 (Nox4) Polyclonal Antibody For Wb, Elisa - Q9Jhi8

Item Number: CM0026775

Introduction

Rabbit polyclonal antibody targeting mouse NOX4, suitable for WB and ELISA. Detects mouse and rat NOX4, a NADPH oxidase generating H₂O₂ and involved in redox signaling and cellular responses.

[Learn More](#)

Parameter	Value
Product Name	NADPH oxidase 4 (NOX4) Rabbit pAb
Remarks/Alias	KOX; KOX-1; RENOX
Species	Mouse
GeneID (Human)	50507
GeneID (Mouse)	50490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 307-424 of mouse NADPH oxidase 4 (NOX4) (NP_056575.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9JHI8
Protein Weight	67 kDa
Shipping	Ice bag