


Anti-Na⁺/K⁺-ATPase Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P05023

Item Number: CM0003232



Introduction

Rabbit monoclonal antibody targeting human Na⁺/K⁺-ATPase alpha-1 subunit (ATP1A1). Validated for WB, IF/ICC, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for ion transport and membrane potential studies.

[Learn More](#)

Parameter	Specification
Product Name	Na ⁺ /K ⁺ -ATPase Rabbit mAb
Remarks/Alias	CMT2DD; HOMGSMR2; Na ⁺ /K ⁺ -ATPase
Species	Human
GeneID (human)	476
GeneID	476
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Na ⁺ /K ⁺ -ATPase (P05023)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05023
Protein Weight	113kDa
Shipping	Ice bag