


Anti-Nanog Rabbit Polyclonal Antibody For Wb, Elisa - Q80Z64

Item Number: CM0019063



Introduction

High-quality mouse NANOG rabbit polyclonal antibody validated for Western blot (WB) and ELISA, with demonstrated human cross-reactivity. Ideal for studying pluripotency and embryonic stem cell self-renewal.

[Learn More](#)

Parameter	Value
Product Name	NANOG Rabbit pAb
Remarks/Alias	ENK; Stm1; ecat4; 2410002E02Rik; NANOG
Species	Mouse
Gene ID (Mouse)	71950
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-305 of mouse NANOG (NP_082292.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q80Z64
Protein Weight	34kDa
Shipping	Ice bag