

Anti-Nap1L1 Monoclonal Antibody For Wb, If/Icc, Elisa - P55209

Item Number: CM0011995



Introduction

Rabbit monoclonal anti-NAP1L1 antibody validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat NAP1L1, a histone chaperone in nucleosome assembly, DNA repair, and neurogenesis.

[Learn More](#)

Parameter	Value
Product Name	NAP1L1 Rabbit mAb
Remarks/Alias	NRP; NAP1; NAP1L; NAP1L1
Species	Human
GeneID	4673
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55209
Protein Weight	45kDa
Shipping	Ice bag