


Anti-Nav1.5/Scn5A Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14524

Item Number: CM0008136



Introduction

Rabbit monoclonal antibody against Nav1.5/SCN5A, validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Targets cardiac voltage-gated sodium channel essential for action potential initiation. Ideal for cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	Nav1.5/SCN5A Rabbit mAb
Remarks/Alias	HB1, HB2, HH1, IVF, VF1, HBBB, ICCD, LQT3, SSS1, CDCD2, CMD1E, CMPD2, PFHB1, Nav1.5
Species	Human
Gene ID (Human)	6331
Gene ID	6331
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 273-357 of human Nav1.5/SCN5A (NP_000326.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14524
Protein Weight	221-227 kDa
Shipping	Ice bag