

Anti-Nbeal2 Polyclonal Antibody For Wb, Elisa - Q6Znj1

Item Number: CM0023962

Introduction

NBEAL2 rabbit polyclonal antibody validated for WB and ELISA in human samples. Ideal for studying neurobeachin-like protein 2 in thrombopoiesis and alpha-granule development. 303kDa, UniProt Q6ZNJ1.

[Learn More](#)

Parameter	Value
Product Name	NBEAL2 Rabbit pAb
Remarks/Alias	GPS; BDPLT4; NBEAL2
Species	Human
GeneID (Human)	23218
GeneID	23218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1800-2100 of human NBEAL2 (NP_055990.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6ZNJ1
Protein Weight	303kDa
Shipping	Ice bag