

Anti-Ncam (Cd56) Rabbit Polymab® Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P13591 / P1595

Item Number: CM0009762



Introduction

Rabbit recombinant monoclonal antibody against NCAM (CD56) for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat; suitable for neuroscience and cancer research.

[Learn More](#)

Parameter	Value
Product Name	NCAM Rabbit PolymAb®
Remarks/Alias	CD56; NCAM; MSK39
Species	Human/Mouse
Gene ID (Human)	4684
Gene ID	4684/17967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-413 of human/mouse CD56 (NP_851996.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13591/P1595
Protein Weight	95kDa
Shipping	Ice bag