


Anti-Ndufa10 (Ci-42Kd) Rabbit Polyclonal Antibody For Wb And Elisa - O95299

Item Number: CM0014236



Introduction

Rabbit polyclonal antibody targeting NDUF10 (CI-42kD), a mitochondrial Complex I accessory subunit. Validated for Western blot and ELISA; reacts with human, mouse, and rat. Suitable for respiratory chain and mitochondrial disease studies.

[Learn More](#)

Parameter	Value
Product Name	NDUF10 Rabbit pAb
Remarks/Alias	CI-42k; CI-42KD; MC1DN22; NDUF10
Species	Human
Gene ID (Human)	4705
Gene ID	4705
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 36-300 of human NDUF10 (NP_004535.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95299
Protein Weight	41kDa
Shipping	Ice bag