


Anti-Ndufb9 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6M9

Item Number: CM0022525



Introduction

Rabbit polyclonal antibody against human NDUF9, an accessory subunit of mitochondrial Complex I. Validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Specification
Product Name	NDUF9 Rabbit pAb
Remarks/Alias	B22; LYRM3; CI-B22; UQOR22; MC1DN24; NDUF9
Species	Human
GeneID (Human)	4715
GeneID	4715
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-170 of human NDUF9 (NP_004996.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6M9
Protein Weight	22kDa
Shipping	Ice bag