

Anti-Neil1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Fi4

Item Number: CM0026092



Introduction

NEIL1 Rabbit pAb (polyclonal antibody) for human, mouse, rat WB and ELISA. Recognizes endonuclease 8-like 1, a DNA repair enzyme involved in base excision repair. Ideal for research on oxidative DNA damage.

[Learn More](#)

Parameter	Value
Product Name	NEIL1 Rabbit pAb
Remarks/Alias	FPG1; NEI1; hFPG1; NEIL1
Species	Human
GeneID (Human)	79661
GeneID	79661
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human NEIL1 (NP_078884.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Fi4
Protein Weight	44kDa
Shipping	Ice bag