


Anti-Nelfe Polyclonal Antibody For Wb, If/Icc, Elisa - P18615

Item Number: CM0029715



Introduction

Rabbit polyclonal antibody against human NELFE, validated for Western blot, immunofluorescence, and ELISA. Suitable for studies on transcriptional pausing, NELF complex dynamics, and HIV-1 latency research.

[Learn More](#)

Parameter	Value
Product Name	NELFE Rabbit pAb
Remarks/Alias	RD; RDP; RDBP; D6S45; NELF-E; NELFE
Species	Human
Gene ID (Human)	7936
Gene ID	7936
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NELFE (NP_002895.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P18615
Protein Weight	43kDa
Shipping	Ice bag