

Anti-Nephrin Polyclonal Antibody For Wb, If, Ihc, Elisa - Q9Qzs7

Item Number: CM0029259

Introduction

Rabbit polyclonal antibody to mouse nephrin (UniProt Q9QZS7). Validated for WB, IF, IHC, ELISA. Reacts with mouse and rat. Suitable for kidney podocyte and glomerular filtration research.

[Learn More](#)

Parameter	Value
Product Name	Nephrin Rabbit pAb
Remarks/Alias	CNF; NPHN; nephrin
Species	Mouse
GeneID (Human)	4868
GeneID	54631
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1033-1256 of mouse Nephrin (NP_062332.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9QZS7
Protein Weight	135kDa
Shipping	Ice bag