

Anti-Neun Mouse Monoclonal Antibody For Wb, If, Ihc, Elisa - A6Nfn3

Item Number: CM0008335

Introduction

High-quality mouse monoclonal antibody targeting NeuN (FOX3/RBFOX3), validated for WB, IF, IHC, ELISA in human, mouse, rat. Splicing regulator protein for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	NeuN Mouse mAb
Remarks/Alias	FOX3; NEUN; FOX-3; HRNBP3
Species	Human
GeneID (Human)	146713
Immunogen	Recombinant protein of human NeuN
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	A6NFN3
Protein Weight	34kDa
Shipping	Ice bag