


Anti-Neurokinin 1 Receptor (Nk1R) Rabbit Monoclonal Antibody For Wb, Elisa - P25103

Item Number: CM0006099



Introduction

Rabbit monoclonal antibody targeting human Neurokinin 1 Receptor (NK1R, TACR1). Suitable for WB and ELISA, cross-reactive with mouse and rat. For tachykinin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Neurokinin 1 Receptor (NK1R) Rabbit mAb
Remarks/Alias	SPR; NK1R; NK1R; TAC1R; Neurokinin 1 Receptor (NK1R)
Species	Human
Gene ID (Human)	6869
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-407 of human NK1R (P25103)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25103
Protein Weight	46 kDa
Shipping	Ice bag