

Anti-Neuropilin-1 (Nrp1) Rabbit Monoclonal Antibody For Wb, Elisa - O14786

Item Number: CM0002082

Introduction

Rabbit monoclonal antibody against human Neuropilin-1 (NRP1, CD304). Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for cardiovascular, angiogenesis, and neuronal research.

[Learn More](#)

Parameter	Value
Product Name	Neuropilin-1 (NRP1) Rabbit mAb
Remarks/Alias	NP1; NRP; BDCA4; CD304; VEGF165R; Neuropilin-1 (NRP1)
Species	Human
GeneID	8829
Immunogen	Synthetic peptide corresponding to amino acids 824-923 of human Neuropilin-1 (NRP1) (O14786)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14786
Protein Weight	103kDa
Shipping	Ice bag