


Anti-Neuropilin-2/Nrp2 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - O60462

Item Number: CM0004959



Introduction

Recombinant rabbit monoclonal antibody targeting human Neuropilin-2 (NRP2). Suitable for Western blot (WB), flow cytometry (FC), and ELISA. Useful for semaphorin/VEGF signaling and cytomegalovirus entry research.

[Learn More](#)

Parameter	Value
Product Name	Neuropilin-2/NRP2 Rabbit mAb
Remarks/Alias	NP2; NPN2; PRO2714; VEGF165R2
Species	Human
Gene ID (Human)	8828
Gene ID	8828
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-855 of human Neuropilin-2/NRP2 (NP_957718.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	O60462
Protein Weight	105kDa
Shipping	Ice bag