

Anti-Nfatc3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q12968

Item Number: CM0027340



Introduction

Rabbit polyclonal antibody against human NFATC3 (NFAT4) for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studying NFAT transcriptional regulation and cardiac development.

[Learn More](#)

Parameter	Value
Product Name	NFATC3 Rabbit pAb
Remarks/Alias	NFAT4; NFATX; NF-AT4c; n339260; NFATC3
Species	Human
Gene ID (Human)	4775
Gene ID	4775
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human NFATC3 (NP_775188.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12968
Protein Weight	116kDa
Shipping	Ice bag