

Anti-Ninj1 Rabbit Polyclonal Antibody For Wb, Elisa - Q92982

Item Number: CM0023781



Introduction

Rabbit polyclonal antibody targeting human Ninjurin-1 (NINJ1). Validated for WB and ELISA; cross-reacts with human, mouse, and rat. Suitable for studies on programmed cell death, inflammation, and nerve injury.

[Learn More](#)

Parameter	Specification
Product Name	NINJ1 Rabbit pAb
Remarks/Alias	NIN1; NINJURIN; NINJ1
Species	Human
Gene ID (Human)	4814
Gene ID	4814
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NINJ1 (NP_004139.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92982
Protein Weight	16kDa
Shipping	Ice bag