


Anti-Nkrf Polyclonal Antibody For Wb, Elisa - O15226

Item Number: CM0027950



Introduction

High-quality rabbit polyclonal antibody against human NKRF (NF-kappa-B-repressing factor) for Western blot and ELISA. Validated in human, mouse, rat. Ideal for transcription regulation research. UniProt O15226.

[Learn More](#)

Parameter	Value
Product Name	NKRF Rabbit pAb
Remarks/Alias	NRF; ITBA4; NKRF
Species	Human
GeneID (Human)	55922
GeneID	55922
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 391-690 of human NKRF (NP_060014.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15226
Protein Weight	78kDa
Shipping	Ice bag