

Anti-Nkx3.1 Rabbit Monoclonal Antibody For Wb And Elisa - Q99801

Item Number: CM0002181



Introduction

High-quality anti-NKX3.1 rabbit monoclonal antibody validated for WB and ELISA. Recognizes human NKX3.1, a prostate-specific tumor suppressor transcription factor. Ideal for prostate cancer research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	NKX3.1 Rabbit mAb
Remarks/Alias	NKX3; BAPX2; NKX3A; NKX3.1
Species	Human
GeneID (human)	4824
GeneID	4824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human NKX3.1 (Q99801)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99801
Protein Weight	26kDa
Shipping	Ice bag