

Anti-Nkx6-1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P78426

Item Number: CM0023661

Introduction

High-quality NKX6-1 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Targets transcription factor essential for beta cell and motor neuron development. Ideal for developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	NKX6-1 Rabbit pAb
Remarks/Alias	NKX6A; NKX6.1; NKX6-1
Species	Human
GeneID (Human)	4825
GeneID	4825
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78426
Protein Weight	38kDa
Shipping	Ice bag