

Anti-Nlrc5 Polyclonal Antibody For Wb, Elisa - Q86Wi3

Item Number: CM0022786

Introduction

Rabbit polyclonal antibody against NLRC5, validated for Western blot and ELISA. Detects human NLRC5 and cross-reacts with mouse. Ideal for studies on NF-kappa-B and interferon signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	NLRC5 Rabbit pAb
Remarks/Alias	NOD4; NOD27; CLR16.1; NLRC5
Species	Human
GenelD (human)	84166
GenelD	84166
Immunogen	Recombinant protein of human NLRC5
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q86WI3
Protein Weight	205kDa
Shipping	Ice bag